

7 Photon Rejuvenation Microcurrent Anti Wrinkle Skin Care Device

Model:LT-007



SAFETY INSTRUCTION

Thank you for purchasing this photon led laser skin care machine

Please read instruction carefully prior to use, and use correctly according to the instruction

Please do not use when you have the following skin problem:

1. Pregnant
2. Using light sensitive medicine to the 560-1200um: Tin AN acids
3. Just after the sun (within 2 weeks)
4. Diabetic
5. Have the sufferer of light allergy
6. Take the recent photosensitive alms
7. The Chronic dissolve a blood suffer

SPECIFICATION

Input voltage: AC100-240V(50/60Hz)

Output Power :15W

Working temperature: -5 centigrade-+40 centigrade

Storage temperature :-5 centigrade-+40 centigrade

PARTS INSTRUCTION



7 Color Photon Therapy Details

Red--One of the most powerful color to increase Blood Circulation. Stimulate skin cell to excrete Collagen to achieve firm & smooth texture.

Green--Anti-aging, Smooth wrinkle & Fine lines, Cure wound

Blue--Excellent effect of Cure acne, Kill bacteria, Help cure sensitive skin

Yellow--Improve skin muscle flexibility, Improve immunity system & neural system

Purple--Relaxing, Improve Lymph Metabolism

Light Blue- Gentle acne treatment, Inflammation skin mitigation

White---Skin whitening, Skin refresh

Features:

During the Photon Beauty treatment process, there is no pain; you will only feel warmth emitted by the pulsing light and a little stimulate form the Microcurrent. This signals that your circulation is being accelerated - those with highly sensitive skin may feel greater stimulation at specific points, but most people will only experience warmth+gentle skin excited.

Operation Steps :

1.BASIC SKIN CARE

STEP 1 Cleanser

Remove Cosmetics and Dead Skin Cells

STEP 2 Lotion

Soften and Moisturize Skin

STEP 3 Photon Beauty+Microcurrent Beauty Device

Activate Skin Cells + Skin Lifting & Firming

STEP 4 Moisturizer

Balance Skin Moisture Content

1. SCARS

STEP 1 Cleanser

Remove Cosmetics and Dead Skin Cells

STEP 2 Lotion

Soften and Moisturize Skin

STEP 3 Photon Beauty+Microcurrent Beauty Device

Activate Skin Cells + Skin Lifting & Firming

STEP 4 Facial Treatment + ULTRASONIC BEAUTY

Provides Nourishment and Stimulates Regeneration of Skin Cells

2.FOLLICULITIS & ACNE

STEP 1 Cleanser

Remove Cosmetics and Dead Skin Cells

STEP 2 Photon Beauty + Microcurrent Beauty Device (Low Level Microcurrent)

Reduce Inflammation +Activate Skin Cells + Skin Lifting & Firming

STEP 3 Medicated Cream

Continue Treatment

3.AGE SPOTS AND PIGMENTATION

STEP 1 Cleanser

Remove Cosmetics and Dead Skin Cells

STEP 2 Lotion

Soften and Moisturize Skin

STEP 3 Photon Beauty+Microcurrent Beauty Device

Activate Skin Cells + Skin Lifting & Firming

STEP 4 Mask

Maintain Moisture Balance

4.FINE LINES AND WRINKLES

STEP 1 Cleanser

Remove Cosmetics and Dead Skin Cells

STEP 2 Lotion

Soften and Moisturize Skin

STEP 3 Photon Beauty+Microcurrent Beauty Device

Activate Skin Cells + Skin Lifting & Firming

STEP 4 Mask + ULTRASONIC BEAUTY

Provides Nourishment and Stimulates Regeneration of Skin Cells

5.SKIN CLARIFICATION

STEP 1 Cleanser

Remove Cosmetics and Dead Skin Cells

STEP 2 Lotion

Soften and Moisturize Skin

STEP 3 Photon Beauty+Microcurrent Beauty Device

Activate Skin Cells + Skin Lifting & Firming

STEP 4 Mask + ULTRASONIC BEAUTY

Provides Nourishment and Stimulates Regeneration of Skin Cells